
KerraContact™ Ag⁺Ag⁺Ag⁺

Advanced silver wound dressing

Information for use

Description

The KerraContact Ag Wound Dressing provides an antibacterial barrier that inhibits bacterial growth in the dressing, as shown *in vitro*, for up to 7 days against *Staphylococcus aureus*, *Enterococcus faecalis*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Acinetobacter baumannii*. The KerraContact Ag Wound Dressing consists of three layers: two non-adherent polyethylene mesh wound contact layers and one absorbent core made of polyester. All three layers are silver coated. The concentration of silver on the KerraContact Ag Wound Dressing is approximately 0.4mg/cm² (2.5% w/w). The KerraContact Ag wound dressing is endotoxin free.

Indications

The KerraContact Ag Wound Dressing is indicated for the management of partial and full thickness wounds, including decubitus ulcers, venous stasis ulcers, diabetic ulcers, first and second degree burns, grafts and donor sites, or other acute or chronic wounds. The dressing may be used over debrided and grafted wounds.

Contraindications

KerraContact Ag wound dressing should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components. Limited information is available regarding the safe use of wound dressings containing silver for pregnant/lactating women, neonates/ small babies and individuals with G6PD deficiency.

Cautions

- Do not use if package is opened or damaged
- Do not re-sterilize. Single use only
- For external use only
- Should not be used on patients with a known sensitivity to silver
- Not compatible with oil-based products, such as petrolatum
- Not compatible with Magnetic Resonance Imaging (MRI) procedures
- Remove dressing prior to administering radiation therapy. A new dressing can be applied following treatment
- Avoid contact with electrode and conductive gels during electronic measurements (e.g. EEG and ECG)
- Store in a dry place (<77°F/25°C). Protect from light
- Avoid exposure to temperatures above 120°F (50°C)
- Permanent skin discoloration may occur after cumulative use of the product at the site of use, where the dressing contacts the skin
- Discontinue use when risk of infection has been lowered. Cumulative use of this product is not intended to exceed 30 days

Directions for Use

- Remove KerraContact Ag Wound Dressing from the package
- Apply the dressing dry or moist and cut to shape as necessary
- Apply the dressing to the wound surface and cover with a secondary dressing as appropriate for the degree of exudate. Secure with surgical tape or a bandage as appropriate
- Examine the dressing daily and change depending on the amount of exudate and the condition of the wound ensuring that a moist environment is maintained
- If the dressing adheres to the wound, moisten or soak the dressing prior to removal
- Avoid forceful removal and disruption of the healing wound
- Dressing maintains inhibition of bacterial growth in the dressing up to 7 days;** as shown *in vitro* against *Staphylococcus aureus*, *Enterococcus faecalis*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Acinetobacter baumannii*

Available Products

KerraContact Ag Wound Dressing is packaged in heat-sealed, peelable, uncoated Tyvek bonded to PerfecFlex opaque white film pouches.

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